


BIOHIT

Innovating for Health

# Unbeatable Accuracy Time after Time



# eLINE

Electronic Single and  
Multichannel Pipettors



# BIOHIT

Innovating for Health

MOVING FROM MANUAL TO  
ELECTRONIC PIPETTING?

# Choose eLINE



### You save time and money

eLINE electronic pipettor is twice as fast as a mechanical pipettor.

### You need only one pipettor instead of two

One Biohit eLINE electronic pipettor covers the volume ranges of two mechanical pipettors.

### Less sick leaves

The electronic operation takes the effort out of pipetting and significantly reduces the risk of Repetitive Strain Injury (RSI).

### You get better results

eLINE electronic pipettor guarantees more accurate and precise results with less effort.



### User independent accuracy

Less experienced users achieve the same precise results as more experienced users.


### Multifunctional tool

eLINE offers a variety of features yet easy to use.

Do you still think that electronic pipetting is expensive? By moving to electronic pipetting you save significantly in the long run!



Comparison of pipetting work with electronic and mechanical pipettors



Both the tip ejection force and the distance thumb has to travel to eject the tips have been measured.

ENJOY YOUR WORK!


## eLINE – the leader in performance & ergonomics

Biohit is the world leader in hand-held electronic pipettors with nearly 20 years of experience in pipettor technologies. Biohit's electronic pipettors have become the standard for electronic pipetting. We put particular emphasis on accuracy and precision, user-friendliness, ergonomics, and efficiency in pipetting.


When you need to pipette and dispense liquids as comfortably as possible, your first choice is eLINE. This award winning pipettor enables you to focus on your work instead of worrying about Repetitive Strain Injury.

### eLINE offers you

- Fully electronic operation modes and tip ejection
- Unbeatable accuracy and precision time after time
- A lot more pipetting modes than in a mechanical pipettor
- 30 times more lighter pipetting than with a mechanical pipettor – most ergonomic pipetting operations on the market
- Applications for example: PCR, ELISA, Volume measurement (manual pipetting mode)



# Double the results!



Multi dispensing  
reduces work stages

*JUST ONE ELINE PIPETTOR BUT*

## A wide range of operation modes


The Biohit eLINE pipettors are available in single-channel models covering the volume range from 0.2  $\mu\text{l}$  to 5000  $\mu\text{l}$  and in 8- and 12-channel models from 0.2  $\mu\text{l}$  to 1200  $\mu\text{l}$ .

- ✗ Pipetting
- ✗ Reverse pipetting
- ✗ Pipetting with mixing
- ✗ Pipetting with cycle counter
- ✗ Manual pipetting with cycle counter\*
- ✗ Multiple dispensing
- ✗ Sequential dispensing
- ✗ Automatic multiple dispensing
- \*Not available in eLINE multi-channel models
- ✗ Diluting
- ✗ Diluting with mixing
- ✗ Multi-aspirating
- ✗ 6 memory settings
- ✗ 9 speeds for aspirating and dispensing
- ✗ Single-point calibration

**“eLINE is the perfect pipettor taking care of both the work and the person using it”**

Lisbeth Lemmetti, Development Analyst, Biomedicine. Blood Centre at Sahlgrenska University Hospital, Sweden

# electronic pipettor



### 1. Easiest tip ejection ever!

- Unique electronic tip ejector function ejects one to 12 tips with the same feather light touch
- Up to 100 times lighter than with a mechanical pipettor
- Tip ejection operation button for both left- and right-handed users

### 2. Perfect tip sealing with Optiload

- Available in multichannel eLINE pipettors
- Allows tip loading with an equal constant force to every channel
- Tips are perfectly sealed to every individual tip cone

### 3. No interruptions because of charging!

- eLINE tells you when the battery is low and you need to recharge it
- You may continue pipetting while charging the battery with AC-adaptor
- When you are not working, keep your eLINE in the charging stand so that the pipettor is always ready to use!
- Charging options: 1-place charging stand, 4-place charging carousel, AC-adaptor



## Unique DC-motor concept for better results

The patented and unique Biohit concept includes the DC-motor, solenoid brake and error control.

### Improved accuracy and precision

The solenoid brake of the eLINE enables over 10 times more precise positioning of the piston compared to a stepper motor driven system. This is a major advantage in multi-dispensing mode where liquid column has to be cut sharply to have even doses. This means superior accuracy and precision.

### Built-in error control for more reliable results

Even the best pipettor may occasionally make errors. Important is to be able to notice them. Biohit electronic pipettors feature a built-in error control. The pipettor controls and monitors the piston movement in real time. If an error occurs while pipetting, the pipettor notifies the user, thus allowing for more reliable results.



## 4. Protect your pipettor and your sample


You can prevent contamination of your pipettor and sample by using replaceable Safe-Cone Filters. These unique filters

- Act as a final barrier to prevent fluids and liquid vapours from reaching the internal components of the pipettor
- Prolong service and lifetime of the pipettor
- Are available for all models > 10 µl
- Safe-Cone Filters must be changed regularly and always in case of over-aspiration

## 10 times more precise dispensing – unbeatable linearity

When comparing the stepper motor technology and the Biohit DC-motor technology for accuracy and precision, especially in the multidispensing mode (d) where even dosing is essential the Biohit DC motor is the best choice. There is no doubt about the superiority of the DC-motor concept.

The DC-motor driven Biohit eLINE and a stepper motor driven pipettor were tested for linearity in the d-mode. The measurement was done 5 times. The DC-motor driven Biohit eLINE allows linear movement, i.e. precise dispensing. The tested stepper motor pipettor generated a linear movement error, which was equal to over  $\pm 1.5 \mu\text{l}$ . This is very alarming as it exceeds the minimum adjustable increment of the pipettor. In the d-mode, where even dosing is important, this kind of an error is significant.




20 x 50 µl linearity measurement in d-mode.



## User-independent accuracy with Biohit eLINE pipettors

Whilst experience and good pipetting practice are essential for reproducible results with mechanical pipettors, the eLINE pipettors produce precise results independent of user experience.

The effect of user experience on pipetting precision



Nine persons were divided into three groups according to their experience in pipetting. Each person pipetted distilled water with the same Biohit mechanical and eLINE pipettors (using P-mode). The imprecision of ten subsequent results of volumes 10, 25, 50, 100, 200 and 1000 µl was recorded and the results are detailed in pictures above.

# eLINE Single and Multichannel Performance specifications and ordering details

| Cat. No. | Channels | Volume Range (µl) | Increment (µl) | Test Volume (µl) | Inacc. (%) | Impr. (%) | Safe-Cone Filters |        | Tip      |
|----------|----------|-------------------|----------------|------------------|------------|-----------|-------------------|--------|----------|
|          |          |                   |                |                  |            |           | Standard          | Plus   |          |
| 73002X   | 1        | 0.2-10 µl         | 0.05           | 10               | 0.90 %     | 0.40 %    | -                 | -      | 10       |
|          |          |                   |                | 5                | 1.00 %     | 0.70 %    |                   |        |          |
|          |          |                   |                | 1                | 2.50 %     | 1.50 %    |                   |        |          |
|          |          |                   |                | 0.2              | 12.00 %    | 10.00 %   |                   |        |          |
| 73004X   | 1        | 5-120 µl          | 0.50           | 120              | 0.40 %     | 0.15 %    | 721008            | 721018 | 300, 350 |
|          |          |                   |                | 60               | 0.60 %     | 0.20 %    |                   |        |          |
|          |          |                   |                | 12               | 2.00 %     | 1.00 %    |                   |        |          |
|          |          |                   |                | 5                | 3.50 %     | 1.50 %    |                   |        |          |
| 73006X   | 1        | 10-300 µl         | 1.00           | 300              | 0.40 %     | 0.15 %    | 721007            | 721017 | 350      |
|          |          |                   |                | 150              | 0.60 %     | 0.20 %    |                   |        |          |
|          |          |                   |                | 30               | 1.50 %     | 0.80 %    |                   |        |          |
|          |          |                   |                | 10               | 3.50 %     | 2.00 %    |                   |        |          |
| 73008X   | 1        | 50-1000 µl        | 5.00           | 1000             | 0.40 %     | 0.15 %    | 721006            | 721016 | 1000     |
|          |          |                   |                | 500              | 0.60 %     | 0.20 %    |                   |        |          |
|          |          |                   |                | 100              | 1.50 %     | 0.50 %    |                   |        |          |
|          |          |                   |                | 50               | 2.00 %     | 1.00 %    |                   |        |          |
| 73010X   | 1        | 100-5000 µl *     | 10.0           | 5000             | 0.50 %     | 0.15 %    | 721006            | 721016 | 5000     |
|          |          |                   |                | 2500             | 0.80 %     | 0.20 %    |                   |        |          |
|          |          |                   |                | 500              | 1.00 %     | 0.40 %    |                   |        |          |
| 73032X   | 8        | 0.2-10 µl         | 0.05           | 10               | 0.90 %     | 0.50 %    | -                 | -      | 10       |
|          |          |                   |                | 5                | 1.50 %     | 0.80 %    |                   |        |          |
|          |          |                   |                | 1                | 4.00 %     | 3.00 %    |                   |        |          |
| 73034X   | 8        | 5-120 µl          | 0.50           | 120              | 0.50 %     | 0.20 %    | 721008            | 721018 | 300, 350 |
|          |          |                   |                | 60               | 0.70 %     | 0.30 %    |                   |        |          |
|          |          |                   |                | 12               | 2.00 %     | 1.50 %    |                   |        |          |
| 73036X   | 8        | 10-300 µl         | 1.00           | 300              | 0.50 %     | 0.20 %    | 721007            | 721017 | 350      |
|          |          |                   |                | 150              | 0.70 %     | 0.30 %    |                   |        |          |
|          |          |                   |                | 30               | 2.00 %     | 1.00 %    |                   |        |          |
| 73039X   | 8        | 50-1200 µl        | 5.00           | 1200             | 0.50 %     | 0.20 %    | 721006            | 721016 | 1200     |
|          |          |                   |                | 600              | 1.00 %     | 0.30 %    |                   |        |          |
|          |          |                   |                | 120              | 2.50 %     | 1.00 %    |                   |        |          |
| 73042X   | 12       | 0.2-10 µl         | 0.05           | 10               | 0.90 %     | 0.50 %    | -                 | -      | 10       |
|          |          |                   |                | 5                | 1.50 %     | 0.80 %    |                   |        |          |
|          |          |                   |                | 1                | 4.00 %     | 3.00 %    |                   |        |          |
| 73044X   | 12       | 5-120 µl          | 0.50           | 120              | 0.50 %     | 0.20 %    | 721008            | 721018 | 300, 350 |
|          |          |                   |                | 60               | 0.70 %     | 0.30 %    |                   |        |          |
|          |          |                   |                | 12               | 2.00 %     | 1.50 %    |                   |        |          |
| 73046X   | 12       | 10-300 µl         | 1.00           | 300              | 0.50 %     | 0.20 %    | 721007            | 721017 | 350      |
|          |          |                   |                | 150              | 0.70 %     | 0.30 %    |                   |        |          |
|          |          |                   |                | 30               | 2.00 %     | 1.00 %    |                   |        |          |
| 73049X   | 12       | 50-1200 µl        | 5.00           | 1200             | 0.50 %     | 0.20 %    | 721006            | 721016 | 1200     |
|          |          |                   |                | 600              | 1.00 %     | 0.30 %    |                   |        |          |
|          |          |                   |                | 120              | 2.50 %     | 1.00 %    |                   |        |          |

X: 0 = Without AC-adaptor 1 = With Universal AC-adaptor (includes Euro, U.S./Jpn, U.K. and China plugs)

\*Note: Minimum volume in P-mode is 500 µl. The 100 µl is possible in d-mode. The specifications are type test specifications achieved under strictly controlled conditions (ISO 8655). Due to Biohit's continuous R&D, specifications may change without prior notice.

## Biohit Tips

| Volume Range (µl) | Packaging    | Qty/Unit     | Cat. No. |
|-------------------|--------------|--------------|----------|
| 0.1-10 µl         | Single Tray  | 10 x 96 Tips | 790010   |
|                   | Refill Tower | 10 x 96 Tips | 790012   |
| 0.5-300 µl        | Single Tray  | 10 x 96 Tips | 790300   |
|                   | Refill Tower | 10 x 96 Tips | 790302   |
| 5-350 µl          | Single Tray  | 10 x 96 Tips | 790350   |
|                   | Refill Tower | 10 x 96 Tips | 790352   |
| 10-1000 µl        | Single Tray  | 10 x 96 Tips | 791000   |
|                   | Refill Pack  | 10 x 96 Tips | 791002   |
| 50-1200 µl        | Single Tray  | 10 x 96 Tips | 791200   |
|                   | Refill Pack  | 10 x 96 Tips | 791202   |
| 100-5000 µl       | Single Tray  | 50 Tips      | 780304   |

All Biohit Single Tray Tips are certified for CE/IVD and RNase, DNase and endotoxin free. More information about tips at [www.biohit.com](http://www.biohit.com).

## Stands & Accessories

| Cat. No. | Item   |
|----------|--|
| 730981   | eLINE Charging Stand for one pipettor, with Universal AC-adaptor (includes Euro, U.S./Jpn, U.K. and China plugs)   |
| 730991   | eLINE Charging Carousel for 4 pipettors, with Universal AC-adaptor (includes Euro, U.S./Jpn, U.K. and China plugs) |
| 725620   | Linear Stand   |
| 731001   | eLINE Replacement Battery  |

Distributed by:



**HEADQUARTERS**  
Biohit Oyj  
Laippatie 1  
00880 Helsinki • Finland  
Tel. +358-9-773 861  
Fax: +358-9-773 86 200  
info@biohit.com  
www.biohit.com

**CHINA**  
Biohit Biotech (Suzhou)  
Co., Ltd.  
Tel: +86-21-6248 5589  
Fax: +86-21-6248 7786  
info.china@biohit.com

**FRANCE**  
Biohit SAS  
Tel: +33-1-3088 4130  
Fax: +33-1-3088 4102  
commercial.france@biohit.com

**GERMANY**  
Biohit Deutschland GmbH  
Tel: +49-6003-82 820  
Fax: +49-6003-828 222  
info@biohit.de

**JAPAN**  
Biohit Japan Co., Ltd.  
Tel: +81-3-5822 0021  
Fax: +81-3-5822 0022  
sales@biohit.co.jp

**RUSSIA**  
Biohit OOO  
Tel: +7-812-327 5327  
Fax: +7-812-327 5323  
main@biohit.ru

**UK**  
Biohit Ltd.  
Tel: +44-1803-315 900  
Fax: +44-1803-315 530  
sales@biohit.co.uk

**USA**  
Biohit Inc.  
Tel: +1-732-922 4900  
Fax: +1-732-922 0557  
pipet@biohit.com

[www.biohit.com](http://www.biohit.com)