

BIOHIT

Innovating for Health

easy
PCR
start



The Biohit Cooling Rack keeps your PCR reagents cool during pipetting. The aluminum rack is autoclavable and it holds 12 micro tubes (1,5 and 2 ml).

Starter Kit Cat. No. 730 820

eLINE 0.2-10 μ l single-channel pipette
eLINE 10-300 μ l single-channel pipette
eLINE Lite electronic dispenser

Charging carousel for 4 pipettes

SafetySpace™ 10 μ l Filter Tip rack, sterile
SafetySpace™ 300 μ l Filter Tip rack, sterile

Dispenser Tip 0.5 ml (10 pcs), sterile
Dispenser Tip 1 ml (10 pcs), sterile

Aluminium cooling rack for 12 micro tubes

Cryo pen (black)

Literature

eLINE® PCR Starter Kit

eLINE PCR Starter Kit offers an easy start in PCR applications. Biohit eLINE electronic pipettes and dispenser speed up your work and guarantee user-independent results and excellent accuracy even in small volumes. Sterile Biohit SafetySpace Filter Tips and Dispenser Tips ensure contamination free pipetting.

eLINE® - unbeatable accuracy

- Fully electronic operation and tip ejection
- User-independent results
- Twice as fast as a mechanical pipette



Autoclavable
lower parts



*"eLINE speeds
up my work!"*

Customer testimonial for eLINE Lite Dispenser

"I do large scale genotyping and DNA/RNA isolation. I dispense the PCR Master Mix to 96 well plate by using the Biohit eLINE Lite Dispenser. The Dispenser makes my work faster and prevents pipetting errors. Also the Biohit Dispenser Tips prevent formation of air bubbles and keep my reagents clean."

*Annariikka Roselli, Engineer of Biotechnology
University of Helsinki, Institute of Biomedicine, Physiology*

www.biohit.com

