



Special Savings on mLINE Starter Kits!

Exclusive Opportunity to Save on mLINE Pipette Purchases



Featuring Biohit's unique Optiload tip mounting mechanism for proper tip sealing

FREE with every purchase!

With every Starter Kit purchased you will receive your choice of one of the following items FREE: 1. mLINE Linear Stand (80087-478) 2. mLINE Carousel Stand (47746-178) or 3. Your choice of two packs of Biohit Pipet Tips. **To select and redeem your free products, go to ca.vwr.com/promotions and search for promo code 3071**



Call your VWR sales representative or 1-800-932-5000 to order. Prices are valid from January 1, 2012 through December 31, 2012.

ORDERING INFORMATION

STARTER KIT: 3 mLINE pipettes

This kit includes 1 single channel pipette for volumes 2-20 µl, 20-200µl & 100-1000µl, a selection of Biohit tips, and a free item of your choice*

VWR Catalog No. MP80078-987

Retail Value: \$1,218.00 Special Price: \$810.84

STARTER KIT: 4 mLINE pipettes

This kit includes 1 single channel pipette for volumes 0.5-10µl, 2-20µl, 20-200µl & 100-1000µl, a selection of Biohit tips, and a free item of your choice*

VWR Catalog No. MP80078-989

Retail Value: \$1,624.00 Special Price: \$1,069.52



- Effective: Jan. 1, 2012 to Dec. 31, 2012
- Free product shipped separately from order, allow 6 to 8 weeks for delivery
- Offer can not be combined with any other promotions, special pricing, or special offers. Multiple redemptions permitted.
- Prices subject to change without notice.
- Complete the online redemption form by January 31, 2013.



TIP REFILL SYSTEM








New Tip Boxes

- The new rigid and strong polypropylene boxes and trays are fully autoclavable.
- The tray holds 96 highest grade non-wettable virgin polypropylene tips in a 8 x 12 micro-plate format.
- The new flip-top boxes are single-hand operated for opening and closing.
- Refill trays snap securely into boxes.
- The color coding of the tip trays allows immediate tip volume recognition.
- Available as non-sterile or pre-sterilized.

The Economic and Ecologic Refill System

- The New Refill System from Biohit is environmentally friendly with completely recyclable packaging.
- The tip trays are easily transferred from a space-saving tower or refill pack to the tip boxes.
- The Refill System range covers the five most popular tip sizes.
- New Value Pack available: 10 refill packs of 9,600, with savings of up to 44%

Biohit ValuePAKs Contain a Total of 9600 Tips

Tip	Vol. Range	Length	Packaging	Type	Qty./Unit	VWR Cat. No.	List Price	Special Price
	0.1 - 10 µl	31.5 mm	Refill Tower	Non Sterile	10x960 Tips 9,600 Tips Total	MP89100-000	\$526.64	\$307.40 <i>(Save \$219.24)</i>
	0.5 - 300 µl	51 mm	Refill Tower	Non Sterile	10x960 Tips 9,600 Tips Total	MP89100-002	\$402.52	\$270.86 <i>(Save \$131.66)</i>
	5 - 350 µl	54 mm	Refill Tower	Non Sterile	10x960 Tips 9,600 Tips Total	MP89100-004	\$462.84	\$273.76 <i>(Save \$189.08)</i>
	10 - 1000 µl	71.5 mm	Refill Pack	Non Sterile	10x960 Tips 9,600 Tips Total	MP89100-006	\$402.52	\$324.80 <i>(Save \$77.72)</i>
	50 - 1200 µl	71.5 mm	Refill Pack	Non Sterile	10x960 Tips 9,600 Tips Total	MP89100-008	\$513.88	\$320.16 <i>(Save \$193.72)</i>

Instructions on How to Use the Refill System

For Refill Tower:



Remove the empty tray from the tip box.



Lift up the upper part of the tower box.



Take one single-tray from the refill tower.



Put the single-tray into the empty tip box.

For Refill Pack:



Detach a Refill Pack.



Open the cover film from the marked corner.



Remove the empty tray from the tip box.



Take the tip tray from the Refill Pack and put it into the empty tip box.



VWR  **1.800.932.5000**
vwr.com

Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in Canadian dollars unless otherwise noted. Offers valid in Canada, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. © 2012 VWR International, LLC. All rights reserved.