


Biohit ValuePak Refill System



**Maximum
Performance,
Maximum
Savings:
BIOHIT
ValuePAK!**

The new Biohit refill system is designed for today's demanding laboratories and scientists. The improved packaging has been designed to be environmentally friendly with reusable tip boxes and completely recyclable tip trays. The result is less plastic waste, greater convenience, and more savings. Tip boxes are tongue-and-groove for easy stacking.

Biohit manufactures and packages all pipet tips using only premium grade materials, following International Organization for Standardization, Declaration of Conformity, and In Vitro Diagnostics standards. All pipet tips are then tested* for fit and function, and packaged completely sealed to ensure that all tips remain clean and RNase/ DNase free until opened. Pipet tips and boxes are fully autoclavable!



BIOHIT



LIMITED TIME OFFER! SAVE UP TO 44%






The Economic and Ecologic Refill System

- The New Refill System from Biohit is environmentally friendly with completely recyclable packaging.
- The tip trays are easily transferred from a space-saving tower or refill pack to the tip boxes.
- The Refill Pack range covers the five most popular tip sizes.
- Tips fit most major pipettors.

New Tip Boxes

- The new rigid and strong polypropylene boxes and trays are fully autoclavable.
- The tray holds 96 highest grade non-wettable virgin polypropylene tips in a 8 x 12 microplate format.
- The new flip-top boxes are single-hand operated for opening and closing.
- Refill trays snap securely into boxes.
- The color coding of the tip trays allows immediate tip volume recognition.

Biohit ValuePAKs Contain a Total of 9600 Tips

Tip	Vol. Range	Length	Packaging	Type	Qty./Unit	VWR Cat. No.	List Price	Special Price
	0.1 - 10 µl	31.5 mm	Refill Tower	Non Sterile	10x960 Tips 9,600 Tips Total	MP89100-000	\$526.64	\$307.40 <i>(Save \$219.24)</i>
	0.5 - 300 µl	51 mm	Refill Tower	Non Sterile	10x960 Tips 9,600 Tips Total	MP89100-002	\$402.52	\$270.86 <i>(Save \$131.66)</i>
	5 - 350 µl	54 mm	Refill Tower	Non Sterile	10x960 Tips 9,600 Tips Total	MP89100-004	\$462.84	\$273.76 <i>(Save \$189.08)</i>
	10 - 1000 µl	71.5 mm	Refill Pack	Non Sterile	10x960 Tips 9,600 Tips Total	MP89100-006	\$402.52	\$324.80 <i>(Save \$77.72)</i>
	50 - 1200 µl	71.5 mm	Refill Pack	Non Sterile	10x960 Tips 9,600 Tips Total	MP89100-008	\$513.88	\$320.16 <i>(Save \$193.72)</i>

Instructions on How to Use the Refill System

For Refill Tower:



For Refill Pack:



- Can not be combined with any other special pricing or offers
- May be cancelled at any time
- Prices subject to change

For information on special pricing, please contact VWR at 1-800-932-5000.

* Tested for chemical resistivity and uniformity & size of orifice.

Call your VWR sales representative or 1-800-932-5000 to order.
Special pricing in effect January 1, 2012 through December 31, 2012.

VWR  1.800.932.5000
vwr.com



Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in Canadian dollars unless otherwise noted. Offers valid in Canada, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. © 2012 VWR International, LLC. All rights reserved.