

BIOHIT

## Limited Time Offer!

Biohit/Pipette MD 4-Step Preventative Maintenance Programs minimize the risk of pipetting errors and prolong the lifespan of your pipette



### 4-Step Mechanical Pipette Program

- Exterior cleaning of pipette
- Replacement of pistons and tip cones
- Leak and Linearity test
- Calibration to manufacturer specifications

PARTS &  
CALIBRATION  
AS LOW AS  
**\$62.00**

Preventative Maintenance Programs are available at highly discounted prices for all Biohit mechanical and electronic pipettes

PARTS &  
CALIBRATION  
AS LOW AS  
**\$99.00**

### 4-Step Electronic Pipette Program

- Exterior cleaning of pipette
- Replacement of battery, pistons, and tip cones
- Leak and Linearity test
- Calibration to manufacturer specifications



Once-a-year replacement of critical parts (tip cones, pistons, and battery for electronic models) with a basic calibration keeps your pipette ahead of the wear and tear curve

PREVENTATIVE MAINTENANCE PROGRAM

# Preventative Maintenance Program

## Programs are available for all Biohit brand pipettes

All pipettes are precision instruments, with mechanical and/or electronic components that are naturally subject to wear and possibly even damage over time, for example:

- The piston surface can deteriorate, wear, or corrode over time from friction and/or chemical attack.
- Tip cones are subject to chemical corrosion and wear, especially in the tip sealing area.
- Batteries within electronic pipettes have a normal lifespan of 1-2 years and must be replaced periodically.



Pipette failure rates can be significantly reduced with the application of a regular Preventative Maintenance Program. Enhance the reliability of your pipette and the integrity of your data!

Sartorius Biohit & Pipette MD recommend Preventative Maintenance programs be applied once a year.

Pipette Type	PM Program Includes	PM Program Cost	Cost if Purchased Separately
1-channel mLINE	Tip cone, piston, calibration	\$62.00	Over a \$130.00 value
8-channel mLINE	Tip cones, pistons, calibration	\$121.00	Over a \$400.00 value
12-channel mLINE	Tip cones, pistons, calibration	\$142.00	Over a \$600.00 value
1-channel eLINE	Tip cone, piston, battery, calibration	\$99.00	Over a \$190.00 value
8-channel eLINE	Tip cones, pistons, battery, calibration	\$157.00	Over a \$575.00 value
12-channel eLINE	Tip cones, pistons, battery, calibration	\$189.00	Over a \$760.00 value

Sartorius Biohit/Pipette MD can provide Preventative Maintenance Programs to any Biohit pipette, including Proline Plus, Proline Electronic, and ePET models. Contact us for additional details.



**For more information:**  
 Biohit Technical Services/Pipette MD  
[www.biohit.com/us](http://www.biohit.com/us)  
[info@pipettemd.com](mailto:info@pipettemd.com)  
 1-877-637-4738

- Offer expires December 31, 2012
- Pricing subject to change without notice
- May not be combined with any other special offers or promotions
- Parts-only pricing without calibration can be quoted separately at a non-promotional level upon request

