



# TRADE-IN EVENT

## mLINE MECHANICAL AND eLINE ELECTRONIC PIPETTES

Biohit and VWR International are now offering an opportunity to trade in ANY BRAND pipette for a Biohit mLINE mechanical or eLINE electronic pipette at a special price. See other side for details, pricing and pipette specifications.

Limited time offer — SAVE OVER 40%



### eLINE FEATURES INCLUDE:

- Direct charging feature allows pipetting while charging. *(Optional one- or four-place charging stands available.)*
- Easy-to-read LCD display informs and instructs the user during pipetting
- Easy-to-use keypad enables fast programming
- Unique electronic one-touch tip ejection
- Electronic start button for all operations
- Attractive ergonomic design means minimal effort
- Volume range color coding for easy pipette selections
- Dispensing head rotates for optimum pipetting convenience
- Individual piston/tip cone assemblies allow easy repair and maintenance
- Optiload tip mounting mechanism for correct tip sealing
- Safe-Cone replaceable tip cone filters help prevent contamination

### mLINE FEATURES INCLUDE:

- 5 color coding caps
- Light plunger button action
- Thumb activated filter ejection
- Ergonomic design and light tip ejection, finger support on the back
- Single, 8- and 12-channel pipettes available
- A clear and easy-to-read display with big numbers and small increments
- A locking mechanism preventing accidental volume changes during pipetting
- Rotating dispensing head for optimal pipetting convenience
- Optiload tip mounting mechanism for correct tip sealing
- Safe-Cone replaceable tip cone filters help prevent contamination



# TRADE-IN EVENT mLINE MECHANICAL AND eLINE ELECTRONIC PIPETTES

Use the special VWR MP catalog numbers to access special pricing on Biohit eLINE and mLINE pipettes.

## eLINE Electronic Pipettes - with NEW universal adapters

Qty:	Pipette w/AC Adapter	Vol. Range	Increment	Inacc. (%)	Impr (%)	Cat No.	List Price	Trade-In Price
_____	eLINE Single Channel	0.1-5 µl	0.05 µl	9.50-1.20	8.00-0.70	<b>MP89202-224</b>	\$887.00	\$564.00
_____	eLINE Single Channel	0.2-10 µl	0.05 µl	0.90-12.00	0.40-10.00	<b>MP89171-638</b>	\$788.26	\$477.00
_____	eLINE Single Channel	5-120 µl	0.50 µl	0.40-3.50	0.15-1.50	<b>MP89171-640</b>	\$788.26	\$477.00
_____	eLINE Single Channel	10-300 µl	1.00 µl	0.40-3.50	0.15-2.00	<b>MP89171-642</b>	\$788.26	\$477.00
_____	eLINE Single Channel	50-1000 µl	5.00 µl	0.40-2.00	0.15-1.00	<b>MP89171-644</b>	\$788.26	\$477.00
_____	eLINE Single Channel	100-5000 µl	10.0 µl	0.50-1.00	0.15-0.40	<b>MP89171-646</b>	\$788.26	\$477.00
_____	eLINE 8-Channel	0.2-10 µl	0.05 µl	0.90-4.00	0.50-3.00	<b>MP89171-648</b>	\$1,382.21	\$836.00
_____	eLINE 8-Channel	5-120 µl	0.50 µl	0.50-2.00	0.20-1.50	<b>MP89171-650</b>	\$1,219.82	\$739.00
_____	eLINE 8-Channel	10-300 µl	1.00 µl	0.50-2.00	0.20-1.00	<b>MP89171-652</b>	\$1,219.82	\$739.00
_____	eLINE 8-Channel	50-1200 µl	5.00 µl	0.50-2.50	0.20-1.00	<b>MP89171-654</b>	\$1,382.21	\$836.00
_____	eLINE 12-Channel	0.2-10 µl	0.05 µl	0.90-4.00	0.50-3.00	<b>MP89171-656</b>	\$1,521.90	\$912.00
_____	eLINE 12-Channel	5-120 µl	0.50 µl	0.50-2.00	0.20-1.50	<b>MP89171-658</b>	\$1,464.65	\$833.00
_____	eLINE 12-Channel	10-300 µl	1.00 µl	0.50-2.00	0.20-1.00	<b>MP89171-660</b>	\$1,370.66	\$833.00
_____	eLINE 12-Channel	50-1200 µl	5.00 µl	0.50-2.50	0.20-1.00	<b>MP89171-662</b>	\$1,513.44	\$912.00

## mLINE Mechanical Pipettes

Qty:	Pipette	Vol. Range	Increment	Inacc. (%)	Impr (%)	Cat No.	List Price	Trade-In Price
_____	mLINE Single Channel	0.1-3 µl	0.002 µl	1.30-10.00%	0.80-6.00%	<b>MP14005-976</b>	\$353.45	\$212.00
_____	mLINE Single Channel	0.5-10 µl	0.01 µl	1.00-5.00%	0.60-4.00%	<b>MP47745-544</b>	\$353.45	\$212.00
_____	mLINE Single Channel	2-20 µl	0.02 µl	0.90-3.00%	0.40-2.00%	<b>MP47745-545</b>	\$353.45	\$212.00
_____	mLINE Single Channel	10-100 µl	0.1 µl	0.80-2.00%	0.15-1.00%	<b>MP47745-546</b>	\$353.45	\$212.00
_____	mLINE Single Channel	20-200 µl	0.2 µl	0.60-2.00%	0.15-0.80%	<b>MP47745-548</b>	\$353.45	\$212.00
_____	mLINE Single Channel	100-1000 µl	1 µl	0.60-1.50%	0.20-0.70%	<b>MP47745-550</b>	\$353.45	\$212.00
_____	mLINE Single Channel	500-5000 µl	10 µl	0.50-0.70%	0.15-0.30%	<b>MP47745-552</b>	\$353.45	\$212.00
_____	mLINE Single Channel	1-10 ml	20 µl	0.60-3.00%	0.20-0.60%	<b>MP89079-110</b>	\$348.46	\$212.00
_____	mLINE 8-Channel	0.5-10 µl	0.01 µl	1.50-6.00%	1.50-6.00%	<b>MP47745-554</b>	\$806.86	\$482.00
_____	mLINE 8-Channel	5-100 µl	0.1 µl	0.70-5.00%	0.25-4.00%	<b>MP47745-556</b>	\$806.86	\$482.00
_____	mLINE 8-Channel	30-300 µl	0.2 µl	0.60-2.00%	0.25-1.00%	<b>MP47745-558</b>	\$806.86	\$482.00
_____	mLINE 12-Channel	0.5-10 µl	0.01 µl	1.50-6.00%	1.50-6.00%	<b>MP47745-560</b>	\$981.71	\$587.00
_____	mLINE 12-Channel	5-100 µl	0.1 µl	0.70-5.00%	0.25-4.00%	<b>MP47745-562</b>	\$981.71	\$587.00
_____	mLINE 12-Channel	30-300 µl	0.2 µl	0.60-2.00%	0.25-1.00%	<b>MP47745-564</b>	\$981.71	\$587.00

First Name \_\_\_\_\_

Last Name \_\_\_\_\_

Job Title \_\_\_\_\_

Division \_\_\_\_\_

Department \_\_\_\_\_

Bldg/Lab/Room: \_\_\_\_\_

Street Address \_\_\_\_\_

P.O. Box/Mail Stop \_\_\_\_\_

City \_\_\_\_\_

State \_\_\_\_\_

Zip Code \_\_\_\_\_

Phone \_\_\_\_\_

Fax \_\_\_\_\_

E-mail \_\_\_\_\_

### Declaration of Decontamination:

I certify that these instruments have been carefully cleaned and decontaminated and pose no danger through bacteriological, virological, chemical, or radioactive contamination. I am aware that I may be held liable for any damages caused by contaminated instruments.

Signature: \_\_\_\_\_

Date: \_\_\_\_\_

Can not be combined with any other special pricing or promotions. May be discontinued without notice. All prices listed in U.S. Dollars. Expires December 31st, 2012. Prices subject to change without notice.



1.800.932.5000 | vwr.com

TIUS2012A

Prices and product details are current when published; subject to change without notice. Certain products may be limited by federal, state, provincial, or local regulations. VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US, void where prohibited by law or company policy, while supplies last. VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owners. Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. © 2012 VWR International, LLC. All rights reserved.