

Receive a **FREE** pipette
with your choice of
SafetySpace™
Filter Tips!
(valued at over \$340)

SafetySpace™ Filter Tips with a free mLINE for ultimate safety

Order any 8 units of SafetySpace™ Filter Tips and get one mLINE single-channel pipette FREE OF CHARGE!



SafetySpace™ Filter Tips for not losing any samples

- Designed with extra space to prevent sample touching the filter
- Can be used safely for foaming liquids, solvents, multiple dispensings and reverse pipetting - for any application that requires extra volume

For premium safety

- E-beam pre-sterilized
- DNase, RNase and endotoxin free certified by lot number
- Manufactured in ISO 8 classified clean room
- CE/IVD certified, ISO 9001 and ISO 13485 compatible

mLINE®

for most ergonomic pipetting

- Prevents injuries and improves pipetting results through very light pipetting forces
- Optiload tip loading mechanism provides perfect tip sealing with low force

HOW TO ORDER

1. Order any 8 units (10x96) of Biohit SafetySpace™ Filter Tips
2. Get one mLINE single-channel mechanical pipette free of charge (ordering details on page 2)

Campaign is valid between November 1, 2011 and February 29, 2012

*To redeem this offer, visit www.vwrsp.com/promotions before 3/31/12, enter Promo Code 3363 in the search box, have your proof of purchase (VWR Invoice or Packing Slip) and fill in the online redemption form.





Ordering details

SafetySpace™ Filter Tips: Choose any 8 units

Tip ID	Volume range	Length	Packaging	Purity Level	Qty/Unit	VWR Cat. No.
L	0.1-10 µl	31.5 mm	Single Tray	Pre-Sterilized	10 x 96	89092-962
N	0.5 - 20 µl	51 mm	Single Tray	Pre-Sterilized	10 x 96	89092-964
O	2 - 120 µl	51 mm	Single Tray	Pre-Sterilized	10 x 96	89092-966
P	5 - 200 µl	54 mm	Single Tray	Pre-Sterilized	10 x 96	89092-968
Q	5 - 300 µl	52.5 mm	Single Tray	Pre-Sterilized	10 x 96	89129-836
S	50 - 1000 µl	78 mm	Single Tray	Pre-Sterilized	10 x 96	89129-838
T	50 - 1200 µl	90 mm	Single Tray	Pre-Sterilized	10 x 96	89129-840

All Single Trays and pre-sterilized Refill Packs are RNase, DNase and endotoxin free


mLINE Single-channel pipettes: Choose any model free of charge

VWR Cat. No.	Channels	Volume Range (µl)	Increment (µl)	Compatible Biohit Tip (µl)	Campaign price
14005-976	1	0.1-3 µl	0.002	10	FREE
47745-544	1	0.5-10 µl	0.01	10	
47745-545	1	2-20 µl	0.02	300	
47745-546	1	10-100 µl	0.10	300, 350	
47745-548	1	20-200 µl	0.20	300, 350	
47745-550	1	100-1000 µl	1.00	1000	
47745-552	1	500-5000 µl	10.0	5000	
89079-110	1	1-10ml	20.0	10000	

See how it works!

Conventional Tip

SafetySpace™ Filter Tip



Air between filter and liquid in reverse pipetting



Prices, product appearance, and specifications are current at the time of printing, subject to change without notice. Availability for certain products may be limited by federal, state, provincial, or local licensing requirements. VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in U.S. dollars unless otherwise noted. Offers valid in USA and Canada, void where prohibited by law or company policy, while supplies last. Visit vwr.com to view our privacy policy and additional disclaimers.

VWR, forms of VWR, and the VWR logo and/or design are either registered trademarks ®, trademarks ™, or service marks ℠ of VWR International, LLC in the United States and/or other countries. All other marks referenced herein are registered trademarks, trademarks, or service marks of their respective owners(s). For a complete list of trademark owners please visit vwr.com.

©2011 VWR International, LLC. All rights reserved.

Campaign is valid between November 1, 2011 and February 29, 2012

1111 Promo Code 3363

